

## **eka-Lpt cushions with cover of incontinence material**

### **CLEANING**

#### **With wipe disinfection:**

To further ensure the product properties of the material we recommend the following products for wipe disinfection:

- Dimozon pur 4%
- Formaldehyd 30%
- Incidin Plus 2%
- Acryl-des 100
- Descogen 10%
- Mikrozid AF liquide 100%
- Perform 6%
- TPH protect 4%

Proceed as follows for wipe disinfection:

The disinfectant used must be used in accordance with the manufacturer's instructions (concentration, exposure time, area of application) to ensure effective disinfection.

**Please observe the manufacturer's instructions!**

For further cleaning instructions, please refer to the label on the product.